

Granulin Rabbit mAb


ABZ3743

Product Overview

Name	Granulin Rabbit mAb
Accession(Primary)	P28799
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW:64 kDa;Observed MW: 55,88 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; IP-1:50-1:100

Target information (P28799)

Gene Name	GRN
About Gene Symbol	GRN (Granulin Precursor) is a cytokine and immune signaling protein involved in immune signaling. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Granulin in 293T lysates
using Granulin antibody.