

CYP27A1 Rabbit mAb


ABZ3745

Product Overview

Name	CYP27A1 Rabbit mAb
Accession(Primary)	Q02318
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:60 kDa;Observed MW: 53kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:50-1:100; FC-1:50-1:100

Target information (Q02318)

Gene Name	CYP27A1
About Gene Symbol	



Western blot analysis of CYP27A1 in mouse liver
lysates using CYP27A1 antibody.