

p73 Rabbit mAb

ABZ3746

Product Overview

Name	p73 Rabbit mAb
Accession(Primary)	O15350
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW:70 kDa;Observed MW: 70,60 kDa
Recommended Dilution	WB-1:500-1:1000; IHC:1:50-1:100; ICC/IF-1:100-1:200; FC-1:50-1:100

Target information (O15350)

Gene Name	TP73
About Gene Symbol	TP73 (Tumor Protein P73) is a transcription factor involved in gene transcription and regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

