

GARNL1 Rabbit mAb


ABZ3747

Product Overview

Name	GARNL1 Rabbit mAb
Accession(Primary)	Q6GYQ0
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:230 kDa;Observed MW: 230 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:50-1:100

Target information (Q6GYQ0)

Gene Name	RALGAPA1
About Gene Symbol	



Western blot analysis of GARNL1 in 293, lysates using GARNL1 antibody.