


RPLP0 Rabbit mAb

ABZ3748

Product Overview

Name	RPLP0 Rabbit mAb
Accession(Primary)	P05388
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:34 kDa;Observed MW: 34 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200

Target information (P05388)

Gene Name RPLP0**About Gene Symbol**

Western blot analysis of RPLP0 in 293, HeLa lysates using RPLP0 antibody.

