

RFC2 Rabbit mAb


ABZ3749

Product Overview

Name	RFC2 Rabbit mAb
Accession(Primary)	P35250
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:39 kDa;Observed MW: 39 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200; FC-1:10-1:100

Target information (P35250)

Gene Name	RFC2
About Gene Symbol	



Western blot analysis of RFC2 in K562, HeLa
lysates using RFC2 antibody.