

PRMT5 Rabbit mAb


ABZ3750

Product Overview

Name	PRMT5 Rabbit mAb
Accession(Primary)	O14744
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:73 kDa;Observed MW: 73 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:20-1:50; FC-1:50-1:100; IP-1:20-1:50

Target information (O14744)

Gene Name	PRMT5
About Gene Symbol	PRMT5 (Protein Arginine Methyltransferase 5) is a methyltransferase involved in chromatin remodeling and epigenetics. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PRMT5 in HepG2, HeLa
lysates using PRMT5 antibody.