

SAMHD1 Rabbit mAb


ABZ3752

Product Overview

Name	SAMHD1 Rabbit mAb
Accession(Primary)	Q9Y3Z3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW:72 kDa;Observed MW: 72 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:10-1:100

Target information (Q9Y3Z3)

Gene Name	SAMHD1
About Gene Symbol	SAMHD1 is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of SAMHD1 in HeLa, MCF7, K562 lysates using SAMHD1 antibody.