

Desmoglein 2/DSG2 Rabbit mAb


ABZ3753

Product Overview

Name	Desmoglein 2/DSG2 Rabbit mAb
Accession(Primary)	Q14126
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:122 kDa;Observed MW: 160 kDa
Recommended Dilution	WB-1:2000-1:5000; IHC-P-1:200-1:500

Target information (Q14126)

Gene Name	DSG2
About Gene Symbol	DSG2 (Desmoglein 2) is a cell junction/adhesion protein involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Desmoglein 2 in HepG2
lysates using Desmoglein 2 antibody.