


HOXA9 Rabbit mAb

ABZ3755

Product Overview

Name	HOXA9 Rabbit mAb
Accession(Primary)	P31269
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW:30 kDa;Observed MW: 30 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (P31269)

Gene Name HOXA9**About Gene Symbol**

Western blot analysis of HOXA9 in U937 lysates
using HOXA9 antibody.

