

Integrin alpha 4 Rabbit mAb


ABZ3758

Product Overview

Name	Integrin alpha 4 Rabbit mAb
Accession(Primary)	P13612
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW:115 kDa;Observed MW: 70,140,150 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P13612)

Gene Name	ITGA4
About Gene Symbol	



Western blot analysis of Integrin alpha 4/CD49D
in Jurkat lysates using Integrin alpha 4/CD49D
antibody.