

Ferritin Heavy Chain Rabbit mAb

ABZ3759


Product Overview

Name	Ferritin Heavy Chain Rabbit mAb
Accession(Primary)	P09528
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW:21 kDa;Observed MW: 21 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:100-1:200

Target information (P09528)

Gene Name Fth1

About Gene Symbol



Western blot analysis of Ferritin Heavy Chain in Mouse liver lysates using Ferritin Heavy Chain antibody.

