

IP3 Receptor Rabbit mAb


ABZ3760

Product Overview

Name	IP3 Receptor Rabbit mAb
Accession(Primary)	Q14643
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:314 kDa;Observed MW: 314 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; IP-1:20-1:50

Target information (Q14643)

Gene Name	ITPR1
About Gene Symbol	



Western blot analysis of IP3 Receptor in Mouse brain lysates using IP3 Receptor antibody.