


Zinc Alpha 2 Glycoprotein Rabbit mAb

ABZ3762

Product Overview

Name	Zinc Alpha 2 Glycoprotein Rabbit mAb
Accession(Primary)	P25311
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW:34 kDa;Observed MW: 41 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IP-1:20

Target information (P25311)

Gene Name AZGP1**About Gene Symbol**

Western blot analysis of Zinc Alpha 2 Glycoprotein in HepG2, Human liver lysates using Zinc Alpha 2 Glycoprotein antibody.

