


Estrogen Related Receptor alpha Rabbit mAb

ABZ3763

Product Overview

Name	Estrogen Related Receptor alpha Rabbit mAb
Accession(Primary)	P11474
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:46 kDa;Observed MW: 50 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:20-1:50; IP-1:20

Target information (P11474)

Gene Name ESRRA**About Gene Symbol**

Western blot analysis of ERR alpha in HeLa
lysates using ERR alpha antibody.

