

Tropomyosin alpha 1 Chain Rabbit mAb


ABZ3764

Product Overview

Name	Tropomyosin alpha 1 Chain Rabbit mAb
Accession(Primary)	P09493
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:33 kDa;Observed MW: 33 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100

Target information (P09493)

Gene Name	TPM1
About Gene Symbol	TPM1 (Tropomyosin 1) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cardiovascular disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Tropomyosin alpha 1 Chain in C6, RAW 264.7 lysates using Tropomyosin alpha 1 Chain antibody.