


CKS2 Rabbit mAb

ABZ3765

Product Overview

Name	CKS2 Rabbit mAb
Accession(Primary)	P33552
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:10 kDa;Observed MW: 10 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100

Target information (P33552)

Gene Name CKS2**About Gene Symbol**

Western blot analysis of CKS2 in Caco-2, 293T lysates using CKS2 antibody.

