

## GIRK-1 Rabbit mAb

ABZ3766

### Product Overview

|                             |                                          |
|-----------------------------|------------------------------------------|
| <b>Name</b>                 | GIRK-1 Rabbit mAb                        |
| <b>Accession(Primary)</b>   | P48549                                   |
| <b>Host Species</b>         | Rabbit                                   |
| <b>Clonality</b>            | Monoclonal                               |
| <b>Isotype</b>              | IgG                                      |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                        |
| <b>Calculated MW</b>        | Calculated MW:57 kDa;Observed MW: 57 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; FC-1:100-1:500         |

### Target information (P48549)

**Gene Name** KCNJ3**About Gene Symbol**

