

Sulfotransferase 2A1 Rabbit mAb


ABZ3767

Product Overview

Name	Sulfotransferase 2A1 Rabbit mAb
Accession(Primary)	Q06520
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:34 kDa;Observed MW: 34 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:500-1:2000; ICC/IF-1:20-1:50; FC-1:50-1:100

Target information (Q06520)

Gene Name	SULT2A1
About Gene Symbol	



Western blot analysis of Sulfotransferase 2A1 in Human fetal liver lysates using Sulfotransferase 2A1 antibody.