


SSRP1 Rabbit mAb

ABZ3768

Product Overview

Name	SSRP1 Rabbit mAb
Accession(Primary)	Q08945
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:81 kDa;Observed MW: 81 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q08945)

Gene Name SSRP1**About Gene Symbol**

Western blot analysis of SSRP1 in K562, 293T
lysates using SSRP1 antibody.

