

AGTR2 Receptor Rabbit mAb

ABZ3769

Product Overview

Name	AGTR2 Receptor Rabbit mAb
Accession(Primary)	P50052
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:41 kDa;Observed MW: 41 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:10-1:100

Target information (P50052)

Gene Name AGTR2

About Gene Symbol



