

TP53INP1 Rabbit mAb


ABZ3770

Product Overview

Name	TP53INP1 Rabbit mAb
Accession(Primary)	Q96A56
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:27 kDa;Observed MW: 27 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200

Target information (Q96A56)

Gene Name	TP53INP1
About Gene Symbol	



Western blot analysis of p53 DINP1 in 3T3 lysates
using p53 DINP1 antibody.