


CD53 Rabbit mAb

ABZ3771

Product Overview

Name	CD53 Rabbit mAb
Accession(Primary)	P19397
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW:24 kDa;Observed MW: 35-55kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; FC-1:50-1:100

Target information (P19397)

Gene Name CD53**About Gene Symbol**

Western blot analysis of CD53 in Raji lysates
using CD53 antibody.

