

S100 alpha 6 Rabbit mAb


ABZ3772

Product Overview

Name	S100 alpha 6 Rabbit mAb
Accession(Primary)	P06703
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW:10 kDa;Observed MW: 10 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:20-1:50; ICC/IF-1:200-1:500

Target information (P06703)

Gene Name	S100A6
About Gene Symbol	



Western blot analysis of S100 alpha 6/PRA in HeLa, A549 lysates using S100 alpha 6/PRA antibody.