

MST3 Rabbit mAb


ABZ3774

Product Overview

Name	MST3 Rabbit mAb
Accession(Primary)	Q9Y6E0
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:49 kDa;Observed MW: 49 kDa
Recommended Dilution	WB-1:500-1:1000; IHC:1:200-1:500; ICC/IF-1:100-1:200; FC-1:20-1:50

Target information (Q9Y6E0)

Gene Name	STK24
About Gene Symbol	



Western blot analysis of MST3 in HeLa lysates
using MST3 antibody.