


FGL1 Rabbit mAb

ABZ3775

Product Overview

Name	FGL1 Rabbit mAb
Accession(Primary)	Q08830
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW:36 kDa;Observed MW: 36 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:10-1:100

Target information (Q08830)

Gene Name FGL1**About Gene Symbol**

Western blot analysis of FGL1 in Human fetal liver lysates using FGL1 antibody.

