

Caspase 14 Rabbit mAb


ABZ3777

Product Overview

Name	Caspase 14 Rabbit mAb
Accession(Primary)	P31944
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:28 kDa;Observed MW: 28 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; IP-1:20; FC-1:50-1:100

Target information (P31944)

Gene Name	CASP14
About Gene Symbol	



Western blot analysis of Caspase 14 in MCF7, HACAT lysates using Caspase 14 antibody.