

COX5A Rabbit mAb


ABZ3778

Product Overview

Name	COX5A Rabbit mAb
Accession(Primary)	P20674
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:17 kDa;Observed MW: 13 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:500-1:2000; IP-1:50-1:100

Target information (P20674)

Gene Name	COX5A
About Gene Symbol	



Western blot analysis of COX5A in HeLa, A431
lysates using COX5A antibody.