


PITX3 Rabbit mAb

ABZ3779

Product Overview

Name	PITX3 Rabbit mAb
Accession(Primary)	O75364
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:32 kDa;Observed MW: 32 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:10-1:100

Target information (O75364)

Gene Name PITX3**About Gene Symbol**

Western blot analysis of PITX3/PTX3 in 293T, U87-MG lysates using PITX3/PTX3 antibody.

