

Cytokeratin 6A Rabbit mAb


ABZ3780

Product Overview

Name	Cytokeratin 6A Rabbit mAb
Accession(Primary)	P02538
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW:60 kDa;Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:20-1:50; ICC/IF-1:20-1:50

Target information (P02538)

Gene Name	KRT6A/KRT6B/KRT6C
About Gene Symbol	



Western blot analysis of Cytokeratin 6 in HACAT
lysates using Cytokeratin 6 antibody.