

HDAC10 Rabbit mAb


ABZ3781

Product Overview

Name	HDAC10 Rabbit mAb
Accession(Primary)	Q969S8
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW:72 kDa;Observed MW: 72 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:100-1:200; FC-1:100-1:200; IP-1:50-1:100

Target information (Q969S8)

Gene Name	HDAC10
About Gene Symbol	HDAC10 (Histone Deacetylase 10) is a histone deacetylase involved in chromatin remodeling and epigenetics. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of HDAC10 in HeLa, 293T, Jurkat, HepG2, 3T3 lysates using HDAC10 antibody.