

## HMGB1 Rabbit mAb


ABZ3782

### Product Overview

|                             |                                                      |
|-----------------------------|------------------------------------------------------|
| <b>Name</b>                 | HMGB1 Rabbit mAb                                     |
| <b>Accession(Primary)</b>   | P09429                                               |
| <b>Host Species</b>         | Rabbit                                               |
| <b>Clonality</b>            | Monoclonal                                           |
| <b>Isotype</b>              | IgG                                                  |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                    |
| <b>Calculated MW</b>        | Calculated MW:25 kDa;Observed MW: 29 kDa             |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-P-1:100-1:200;<br>FC-1:20-1:100 |

### Target information (P09429)

|                          |       |
|--------------------------|-------|
| <b>Gene Name</b>         | HMGB1 |
| <b>About Gene Symbol</b> |       |



Western blot analysis of HMGB1 in HeLa lysates  
using HMGB1 antibody.