

Flightless I Rabbit mAb


ABZ3785

Product Overview

Name	Flightless I Rabbit mAb
Accession(Primary)	Q13045
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:145 kDa;Observed MW: 145 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200

Target information (Q13045)

Gene Name	FLII
About Gene Symbol	



Western blot analysis of Flightless I in HeLa, A549, A431 lysates using Flightless I antibody.