

Dishevelled 3 Rabbit mAb


ABZ3786

Product Overview

Name	Dishevelled 3 Rabbit mAb
Accession(Primary)	Q92997
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:78 kDa;Observed MW: 95 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:20-1:50; FC-1:100-1:200

Target information (Q92997)

Gene Name	DVL3
About Gene Symbol	



Western blot analysis of Dishevelled 3 in K562
lysates using Dishevelled 3 antibody.