

## Ly6g Rabbit mAb

ABZ3788

### Product Overview

|                             |                                    |
|-----------------------------|------------------------------------|
| <b>Name</b>                 | Ly6g Rabbit mAb                    |
| <b>Accession(Primary)</b>   | P35461                             |
| <b>Host Species</b>         | Rabbit                             |
| <b>Clonality</b>            | Monoclonal                         |
| <b>Isotype</b>              | IgG                                |
| <b>Species Reactivity</b>   | Mouse                              |
| <b>Calculated MW</b>        | Calculated MW: 14kDa               |
| <b>Recommended Dilution</b> | IHC-P-1:200-1:2000; FC-1:200-1:500 |

### Target information (P35461)

|                          |      |
|--------------------------|------|
| <b>Gene Name</b>         | Ly6g |
| <b>About Gene Symbol</b> |      |