

## Cleaved C-terminal GSDMD Rabbit mAb

ABZ3790

### Product Overview

|                             |                                     |
|-----------------------------|-------------------------------------|
| <b>Name</b>                 | Cleaved C-terminal GSDMD Rabbit mAb |
| <b>Accession(Primary)</b>   | Q9D8T2                              |
| <b>Host Species</b>         | Rabbit                              |
| <b>Clonality</b>            | Monoclonal                          |
| <b>Isotype</b>              | IgG                                 |
| <b>Species Reactivity</b>   | Mouse                               |
| <b>Calculated MW</b>        | Calculated MW: 22 kDa               |
| <b>Recommended Dilution</b> | WB-1:500-1:1000                     |

### Target information (Q9D8T2)

|                          |       |
|--------------------------|-------|
| <b>Gene Name</b>         | GSDMD |
| <b>About Gene Symbol</b> |       |