

GSDMD Rabbit mAb

ABZ3792

Product Overview

Name	GSDMD Rabbit mAb
Accession(Primary)	Q9D8T2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 53 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:200-1:1000; FC-1:200-1:500; IP-1:20-1:50

Target information (Q9D8T2)

Gene Name	GSDMD
About Gene Symbol	