

PI3 Kinase p110 delta Rabbit mAb

ABZ3793


Product Overview

Name	PI3 Kinase p110 delta Rabbit mAb
Accession(Primary)	O00329
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 119 kDa; Observed MW: 119 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200

Target information (O00329)


Gene Name PIK3CD

About Gene Symbol



Western blot analysis of PI3-Kinase p110 delta in Jurkat, Raji, K562, RAW 264.7, Rat brain, Rat kidney lysates using PI3-Kinase p110 delta antibody.

Copyright© 2024 Zen-bioscience Limited.



Western blot analysis of PI 3 Kinase p110 delta in (1) K562 lysates; (2) RAW 264. lysates; (3) rat kidney lysates using PI3 Kinase p110 delta antibody.