


ApoH Rabbit mAb

ABZ3794

Product Overview


Name	ApoH Rabbit mAb
Accession(Primary)	P02749
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 38 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:100

Target information (P02749)

Gene Name APOH**About Gene Symbol**

Western blot analysis of Apolipoprotein H in Human embryo lysates using Apolipoprotein H antibody.

Copyright© 2024 Zen-bioscience Limited.



Western blot analysis of ApoH in Human plasma lysates using ApoH antibody.