

FGFR2 Rabbit mAb

ABZ3797


Product Overview

Name	FGFR2 Rabbit mAb
Accession(Primary)	P21802
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 92 kDa; Observed MW: 145 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:10-1:100

Target information (P21802)


Gene Name FGFR2

About Gene Symbol



Western blot analysis of FGFR2 in Jurkat□MCF7 lysates using FGFR2 antibody.

Copyright© 2024 Zen-bioscience Limited.



Western blot analysis of FGFR2 in (1) MCF-7 lysates; (2) mouse brain lysates using FGFR2 antibody.