


CPT2 Rabbit mAb

ABZ3798


Product Overview

Name	CPT2 Rabbit mAb
Accession(Primary)	P23786
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 74 kDa; Observed MW: 74 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200


Target information (P23786)

Gene Name CPT2**About Gene Symbol**


Copyright© 2024 Zen-bioscience Limited.



Western blot analysis of CPT2 in MCF-7 lysates using CPT2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human ovarian cancer using CPT2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of CPT2 in MCF-7 using CPT2 antibody.

