

HAPLN1 Rabbit mAb


ABZ3799

Product Overview

Name	HAPLN1 Rabbit mAb
Accession(Primary)	P10915
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 40 kDa; Observed MW: 48 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:20-1:50; ICC/IF-1:100-1:200; FC-1:200-1:2000


Target information (P10915)

Gene Name	HAPLN1
About Gene Symbol	




Western blot analysis of HAPLN1 in MCF7,Caco2 lysates using HAPLN1 antibody.

Copyright© 2024 Zen-bioscience Limited.



Western blot analysis of HAPLN1 in (1) Caco2 lysates; (2) mouse spleen lysates using HAPLN1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human kidney using HAPLN1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.