

TOP2A Rabbit mAb


ABZ3800

Product Overview

Name	TOP2A Rabbit mAb
Accession(Primary)	P11388
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 174 kDa; Observed MW: 174 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IP-1:50


Target information (P11388)

Gene Name	TOP2A
About Gene Symbol	TOP2A (DNA Topoisomerase II Alpha 2) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of Topoisomerase II alpha in Jurkat lysates using Topoisomerase II alpha antibody.



Western blot analysis of Topoisomerase II alpha in (1) HeLa lysates; (2) Jurkat lysates using Topoisomerase II alpha antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Topoisomerase II alpha antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.