

CD168 Rabbit mAb


ABZ3803

Product Overview

Name	CD168 Rabbit mAb
Accession(Primary)	O75330
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 84 kDa; Observed MW: 84 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; IP-1:10-1:100


Target information (O75330)

Gene Name	HMMR
About Gene Symbol	HMMR (Hyaluronan Mediated Motility Receptor) is a membrane receptor involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of CD168 in T47-D lysates using CD168 antibody.



Western blot analysis of CD168 in LnCaP lysates using CD168 antibody.



Immunohistochemistry analysis of paraffin-embedded Human testis using CD168 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.