

## CD21 Rabbit mAb


ABZ3804

### Product Overview

|                             |                                                                           |
|-----------------------------|---------------------------------------------------------------------------|
| <b>Name</b>                 | CD21 Rabbit mAb                                                           |
| <b>Accession(Primary)</b>   | P20023                                                                    |
| <b>Host Species</b>         | Rabbit                                                                    |
| <b>Clonality</b>            | Monoclonal                                                                |
| <b>Isotype</b>              | IgG                                                                       |
| <b>Species Reactivity</b>   | Human                                                                     |
| <b>Calculated MW</b>        | Calculated MW: 113 kDa; Observed MW: 155 kDa                              |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:1000-1:5000;<br>ICC/IF-1:20-1:50; IP-1:50-1:100 |


### Target information (P20023)

|                          |     |
|--------------------------|-----|
| <b>Gene Name</b>         | CR2 |
| <b>About Gene Symbol</b> |     |




Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of CD21 in Daudi,Raji  
lysates using CD21 antibody.



Western blot analysis of CD21 in Raji lysates  
using CD21 antibody.



Immunohistochemistry analysis of paraffin-embedded Human spleen using CD21 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.