

IRF2 Rabbit mAb


ABZ3805

Product Overview

Name	IRF2 Rabbit mAb
Accession(Primary)	P14316
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 39 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:20-1:50; ICC/IF-1:200-1:500


Target information (P14316)

Gene Name	IRF2
About Gene Symbol	



Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of IRF2 in HeLa and NIH/3T3 lysates using IRF2 antibody.



Western blot analysis of IRF2 in (1) HeLa lysates; (2) 3T3 lysates using IRF2 antibody.