

ATF1 Rabbit mAb

ABZ3806


Product Overview

Name	ATF1 Rabbit mAb
Accession(Primary)	P18846
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 29 kDa; Observed MW: 36 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:10-1:100

Target information (P18846)


Gene Name ATF1

About Gene Symbol



Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of ATF1 in K562, HeLa, HepG2 lysates using ATF1 antibody.



Western blot analysis of ATF1 in Hela lysates
using ATF1 antibody.