

Urokinase Rabbit mAb

ABZ3807


Product Overview

Name	Urokinase Rabbit mAb
Accession(Primary)	P00749
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 49 kDa; Observed MW: 49 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:100-1:200

Target information (P00749)


Gene Name PLAU

About Gene Symbol



Western blot analysis of uPA in PC-3, LnCaP lysates using uPA antibody.

Copyright© 2024 Zen-bioscience Limited.



Western blot analysis of Urokinase in LnCaP
lysates using uPA antibody.