

HLA-DPB1 Rabbit mAb


ABZ3808

Product Overview

Name	HLA-DPB1 Rabbit mAb
Accession(Primary)	P04440
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 29 kDa; Observed MW: 29 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:2000; ICC/IF-1:50-1:200; IP-1:10-1:100


Target information (P04440)

Gene Name	HLA-DPB1
About Gene Symbol	



Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of MHC Class II in Human tonsil lysates using MHC Class II antibody.



Western blot analysis of MHC Class II in Daudi lysates using HLA DPB1 antibody.