

Claudin 4 Rabbit mAb


ABZ3811

Product Overview

Name	Claudin 4 Rabbit mAb
Accession(Primary)	O14493
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 22 kDa; Observed MW: 18 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:500-1:2000; IP-1:20-1:50


Target information (O14493)

Gene Name	CLDN4
About Gene Symbol	



Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of Claudin 4 in PC-3,SW480,MCF-7 lysates using Claudin 4 antibody.



Western blot analysis of Claudin 4 in MCF-7 lysates using Claudin 4 antibody.